BRAUN

Oral-B

sonic complete

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Safety Precautions

- Please read the «Safety Precautions» thoroughly beforehand to ensure proper use of the product.
- The below precautions are to ensure that you use the product safely and correctly, so as to prevent injury or harm to yourself or others. This section has been divided into two categories «Warnings» and «Precautions» to clearly indicate the severity of the harm or damage that may potentially arise from misuse of this product.

Each precaution contains important information relating to your safety, so please make sure you follow them.

Warning: Pertains to information relating to potential death or severe injury

Precaution: Pertains to information relating to potential injury or physical harm

Examples of Graphic Indicators

- △ Symbol highlights Warning or Precautions. Specific precautions are pictured inside the symbol. For example, the symbol on the left represents precaution to be taken to avoid electrocution.

- ⊘ Symbol indicates a forbidden action. Forbidden actions are pictured inside or near the symbol. For example, the symbol on the left alerts that disassembly is forbidden.


- Symbol indicates certain recommended actions. Specific information is pictured inside the symbol. For example, the symbol on the left indicates that the product should be unplugged from the power socket in a certain manner.

Warnings

- Individuals other than certified repair technicians should not disassemble or attempt to repair the product. Please do not attempt to modify the product yourself. This could lead to fire, electrocution, or injury. For repairs, please speak with a representative from the purchase location or with the Customer Service Center.

- Please assure that the product is unplugged during product maintenance. In addition, please do not insert or remove the plug from the power socket with wet hands, since this could lead to electrocution or injury.

Precautions

- When removing the plug from the power socket, please pull the plug and avoid pulling directly on the power cord. This will protect against electrocution or a short circuit, which can lead to the outbreak of fire.

- Please check that the toothbrush is not shaky or damaged before use, since this could harm your gums or lips.

- Please do not use the product if the power cord or plug is damaged or if the plug does not securely fit into the power socket. This will protect against electrocution or a short circuit, which can lead to the outbreak of fire.

- Please do not damage, break, convert, forcibly bend, stretch, twist or bundle together the power cord. In addition, placing a heavy object on the cord, lodging the cord between objects, or converting the power cord could cause damage to the cord. A damaged cord could lead to fire or electrocution.


- Please do not use with a current outside the range AC 100-120V. Such action could lead to fire or electrocution.

- Please do not press firmly on the brush head. This could damage your gums or lips.

- Please do not immerse or pour water on the charging station. Such action could lead to a short circuit or electrocution.


After reviewing the above, please assure that this document is easily accessible to those using the product.

Index

Safety Precautions	14
Braun Oral-B Sonic Complete™ Specifications	17
The Features of the Braun Oral-B Sonic Complete	18
Charging the Braun Oral-B Sonic Complete	19
The Recommended Brushing Technique.....	20
Customize your Clean.....	21
Caring for the Braun Oral-B Sonic Complete	22
Changing the Brush Head.....	22
Caring for the Environment.....	23
Braun Oral-B Sonic Complete Q&A.....	24
Guarantee.....	24
中文	25

Braun Oral-B Sonic Complete Specifications


Every **Braun Oral-B Sonic Complete** features the following:

- 1) Brush head
- 2) ON/OFF switch
- 3) Sonic handle
- 4) Brushing modes
- 5) Charge indicator
- 6) Charger base
- 7) Storage tower for additional brush heads
- 8) Wall holder

Motion	Frequency
	Clean 31,000 times/minute
	Soft 29,000 times/minute
	Massage 24,000-33,600 times/minute
Power supply	Rechargeable
Recharge time	16 hours
Rating	AC 100-120V, 50-60Hz, 2W
Usage time	2 weeks (approximately 56 minutes)
Measurements	Main unit: diameter 36mm x height 196mm (including brush head, 243mm) Charger base: width 47mm x depth 110mm x height 178mm
Mass	143g (main unit) 238g (charger base including power cord)
Accessories	Brush head, wall holder, screws, travel case*, travel cap*

*S18 535 3 and S1800 only

How to attach the wall holder



- Use the enclosed screws to attach fixture to wall.
- Place extra length of power cord in the power cord storage space in the back of the charger station. Slide charger base into place from above and secure.
- If screws cannot be secured into the wall, use plastic screws before inserting screws.

The Features of the Braun Oral-B Sonic Complete

Braun Oral-B Sonic Complete has the following features designed to offer you complete mouth care:

Premium bristling technology

- 1) CrissCross™ bristles
Advanced multi-textured bristles proven to better clean tooth surfaces, between teeth and along the gumline.
- 2) Interdental bristles
Raised for better interdental cleaning.
- 3) Power Tip® bristles
Extra-long bristles at the tip to help clean hard-to-reach areas.
- 4) Indicator® bristles
Blue color fades halfway with use, signaling the need to replace the brush head.


Microsonic brushing modes for complete mouth care

- «Clean»
for a complete cleaning of teeth and gums.
- «Soft»
for the gentle care of the tongue and other sensitive areas.
- «Massage» (S18 535 3 and S1800 only)
for a gentle stimulation designed to help keep gums firm and healthy.

Timers to encourage thorough brushing

- **Professional Timer:** A short stuttering sound at 30-second intervals reminds you to brush the four quadrants of your mouth equally for a better overall clean.
- **2-Minute Timer:** A longer stuttering sound indicates when the recommended 2-minute brushing time has been reached.

Charging the Braun Oral-B Sonic Complete

- When you place your handle on the charger base and plug the base into an electrical outlet, the blue charge indicator will light up. While charging, the light blinks. Once the handle is fully charged, the light will stop blinking and shine continuously. If you then use your brush and place it back on the charger, the charging light will again blink to let you know it is charging.
- A full charge lasts for approximately two weeks of regular brushing, so you can easily take your **Braun Oral-B Sonic Complete** with you on vacation.
- Your **Braun Oral-B Sonic Complete** is designed to provide you with an effective clean consistently. If the battery has run too low to provide this clean, your brush will stop.
- We recommend keeping your brush on the charger base when not in use to maintain a full battery charge. If your rechargeable battery does run down completely, it will take about 16 hours for your **Braun Oral-B Sonic Complete** to fully charge.

When removing the plug from the power socket, please make sure that you pull the plug and not just the power cord. Pulling just the cord could result in electrocution or a short circuit, which could lead to the outbreak of fire.

- It is recommended that you (1) Charge and (2) Discharge once every 6 months.
 - (1) Fully charge for 16 hours, then use until completely discharged.
 - (2) When the unit is completely discharged, return unit to charger base and fully charge for 16 hours.

Please do not use below or above AC 100-120V. This could lead to fire or electrocution.

- Because the charger base uses electromagnetic induction, please do not place on top of stainless steel or other metallic surfaces when in use. Recharge capacity may be affected and recharging may become difficult.
- Please check that the power cord of the charger base is not damaged. If there are signs of damage, please refrain from use and contact the purchase location or the local Braun Repair Center.

- Rechargeable batteries may grow weaker over time. A notable reduction in running time, even after charging, may mean that the battery has begun to deteriorate. If repairs are required (charges apply), please see your nearest sales outlet or Braun Repair Center.

Please do not damage, break, convert, forcibly bend, stretch, twist or bundle together the power cord.

Individuals other than certified repair technicians should not disassemble or attempt to repair the product. Please do not attempt to modify the product yourself. This could lead to fire, electrocution, or injury.

The Recommended Brushing Technique



1. Attach brush head firmly to the main unit.
 - A 2mm gap exists between the brush head and the main unit. This space is necessary for the movement of the brush head. This is not a product defect.

Please check that the toothbrush is not shaky or damaged before use, since this could harm your gums or lips.

The brush head may break or become damaged if a strong force is exerted on it. Please refrain from using brush head in this situation.

2. Please turn switch ON after placing the toothbrush in your mouth.
 - When toothpaste is used, please be particularly careful when switching product ON to prevent the toothpaste from spraying.
 - Angle the brush along the gumline.
 - Hold the handle lightly in your hands and turn switch ON.
3. Learn how to brush effectively.
 - Gently move the brush head in a slightly circular motion. Allow the longer interdental bristles to reach between your teeth, while the CrissCross bristles focus on cleaning the tooth surface and gumline regions (see diagram).

- After a few seconds, glide the bristles to the next section and pause while you gently circulate the bristles again. Try using only light pressure.
- Adjust your brushing, so that you keep as many bristles in contact with your teeth as possible. Don't forget to reach all areas, including the inner and chewing surfaces, and behind the back teeth. Remember to floss daily.

There is no need to press firmly on the brush head, since such action could damage your gums or lips.

4. Toothpaste

This product can be used with the toothpaste of your choice. Please wet brush before applying toothpaste.

Customize your Clean

Microsonic brushing modes for complete mouth care:

- «Clean» for a complete cleaning of teeth and gums.
- «Soft» for the gentle care of the tongue and other sensitive areas.
- «Massage» (S18 535 3 and S1800 only) for a gentle stimulation designed to help keep gums firm and healthy.

How to switch between modes:

- To start your **Braun Oral-B Sonic Complete**, push the on/off button once, and your toothbrush automatically begins vibrating in the «Clean» mode.
- To shift from the «Clean» mode to the «Soft» mode, push the on/off button once more.
- If your **Braun Oral-B Sonic Complete** has a «Massage» mode (S18 535 3 and S1800), pushing the on/off button again will bring you from «Soft» to «Massage» mode.
- To turn your brush off from the «Massage» mode, push the on/off button again.
- If you wish to turn your **Braun Oral-B Sonic Complete** off in the «Clean» or «Soft» mode, push and hold down the on/off button.

You may experience a slight tickling or tingling sensation when you first use **Braun Oral-B Sonic Complete**. This will subside as your teeth and gums become accustomed to the **Braun Oral-B Sonic Complete** cleaning experience. If you experience this initial

- Clean
- Soft
- Massage


tingling, you may want to brush in the «Soft» mode at first.

During the first few days of using your **Braun Oral-B Sonic Complete**, it is possible that your gums may bleed slightly. This should stop after a few days. However, should it persist, please consult your dentist.

Caring for the Braun Oral-B Sonic Complete

Constant moisture or toothpaste on the handle or brush head may lead to moldiness, so please follow these simple cleaning steps daily.

- Turn switch OFF, remove brush head from handle, and rinse in running water as shown in the diagram on the left.
- Please use a cloth to wipe off any moisture or dust from the handle or charger base. Please clean regularly.
- Please remove brush from handle and dry in well-ventilated place.
- The brush head storage area in the charger base is removable. Please rinse under running water.
- When placing the storage area back into the charger base, open lid to the upright position, and insert storage area from the above by following the markers.


Please do not immerse or pour water onto charger base, since there could be a danger of a short circuit or electrocution.

Changing the Brush Head

The blue Indicator bristles fade halfway to signal when it is time to replace the brush head. It is recommended that brush heads be replaced every 3 to 4 months. A worn-out brush head is not only less effective at removing plaque than a new brush head, it can also damage your gums. When brushing twice a day for two minutes each time, the blue Indicator bristles will fade halfway in approximately 3 to 4 months, signaling that it is time to replace the brush head.

* Replacement brush heads are sold separately. Please purchase at your nearby store.

Please note that brush heads other than those designed for **Braun Oral-B Sonic Complete** (SR18-2) should not be used.

Caring for the Environment


- This symbol indicates that the **Braun Oral-B Sonic Complete** contains a nickel cadmium battery. At the end of the product's life, please dispose of the rechargeable battery inside the handle properly. This battery cannot be replaced, but it is easily removed for recycling or proper disposal.
- To help protect your environment, please dispose of the battery properly at approved recycling or disposing locations in accordance with local regulations.

For Recycling Locations

For details, please see the Japan Battery Recycling Center homepage.

Homepage: <http://www.JBRC.com>

Please do not disassemble the handle except for disposal purposes. Once disassembled, repair of the handle is not possible.

- ① Bring bottom of handle in line with the groove on the bottom of the charger base and turn counter-clockwise for approximately 60 degrees (until it can no longer be turned).
- ② Point handle downwards and push against a hard surface. The motor and internal parts containing the batteries will come out from the handle.
- ③ Pull out internal parts.
- ④ Being careful not to injure yourself, take out the motor and cylindrical batteries.
- ⑤ Remove batteries by using a screwdriver or the like to break off the end of the batteries from the printed circuit board.
- ⑥ Please take the used batteries to the recycling box of your nearest recycling location.

Braun Oral-B Sonic Complete Q&A

My teeth and gums are sensitive.

How should I best brush with the Braun Oral-B Sonic Complete?

We recommend that people with sensitive teeth or gums start by brushing in the «Soft» mode for 1 to 2 weeks.

What should I do to effectively remove coffee, tea, and tobacco stains from my teeth?

Braun Oral-B Sonic Complete is proven to remove stains for naturally whiter teeth. For heavy stains, we recommend brushing slightly longer to aid in their removal.

My dentist told me that I must take better care of my gums.

Is Braun Oral-B Sonic Complete effective?

Yes. Clinical trials have shown that

Braun Oral-B Sonic Complete can help maintain healthy, firm gums. Proper brushing for the right amount of time is key for this to occur. **Braun Oral-B Sonic Complete** includes both the Professional Timer as well as the 2-Minute Timer. The Professional Timer lets you know when to move to the next quadrant of your mouth, so that you brush all four quadrants of your mouth equally and the 2-Minute Timer signals when the recommended brushing time has passed.

Can Braun Oral-B Sonic Complete be used by children?

Yes. Children should use the product under parental supervision to ensure that all teeth and areas of the mouth are properly cleaned. We recommend using the «Soft» mode of brushing on children's teeth, as this mode is specially designed for more delicate tooth and gum areas.

Someone in my family wears braces. What is the recommended way to brush with Braun Oral-B Sonic Complete?

Gently place the brush at an angle between the brackets and gums. Using the «Clean» mode, move the brush in a light circular motion to clean both the brackets and teeth.

I have bridges, crowns and fillings in my mouth. What is the best way to brush with my Braun Oral-B Sonic Complete?

Like all Braun Oral-B power toothbrushes, **Braun Oral-B Sonic Complete** can be used by people with bridges, implants and crowns. We recommend using the «Clean» mode to effectively remove plaque and bacteria for healthy teeth and gums.

What should I do if I have comments or questions about the Braun Oral-B Sonic Complete?

Please send your comments and questions to Braun's Customer Services Center (Japan tel. 0120-13-6343, Taiwan tel. 0800-005-998). You can also contact us through our homepage (www.oralb.com or www.braun.co.jp).

Guarantee

We grant 1 year guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose. This guarantee extends to every country where this appliance is supplied by Braun or its appointed distributor.

This guarantee does not cover: damage due to improper use, normal wear or use, especially regarding the brush heads, as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Braun Repair Center.

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